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YC-019 pilot-scale spray dryer is an innovative product of Shanghai Pilotech Instrument & Equipment Co., Ltd. through years of exploration and research. Compared to current pilot-scale spray dryers, this machine has the advantages of small volume, lower noise, low power, high recovery rate and convenient operation, which can be used in the laboratory. The whole machine occupies less than one square meter. It is configured with one-button start-up and the color LCD touch screen for easy operation. Also, it can run automatically or manually for monitoring in the experimental process. It generates spray particles with a diameter of 1 to 100 microns and performs excellent in drying.

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Two modes of run: Automatic-mode, and Eye-monitored mode for the purpose of easily controlling experimental process .

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The spray head is a concentric one. It guarantees no eccentricity during atomization that may lead to spraying to the bottle wall. The spray head can move up and down after installed to facilitate the adjustment of the atomization position to improve the spray drying

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To assure accurate temperature control and easy parameter adjusting, YC-019 applies PID. Temp. Precision: $\pm 1^{\circ}\text{C}$, dry cycle time: 1.0-2.5S

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Temperature protection: At the end of the experiment, the heater temperature is extremely high, which requires the blower continue working to reduce the internal temperature to ensure the safety of the equipment. YC-019 Laboratory Spray Dryer's program can automatically control the rotation of the blower. Even if the operator makes a mistake, user still cannot shut the blower down until the system temperature is reduced to the system default safe state. Then the blower will automatically stop;

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As standard a touch screen to display temperature, airflow volume, air pressure, pump speed and de-blocker frequency.

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Standard with 2.0 mm jet and other sizes are available as accessories

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